

USER MANUAL



AIRPORT EMERGENCY TRANSPORT TROLLEY

For model
B6-030-4000

NOTE!

**The operator must read the user manual
carefully BEFORE the product is used,
to avoid unnecessary risks to the patient and operator.**

TABLE OF CONTENTS

Foreword

Warranty

Important information

1. Description of functions

1.1 General description

2. Operation

2.1 Height adjustment

2.2 The castors

2.3 The slide

2.4 Towards the airplane

3. Maintenance & Cleaning

3.1 Maintenance, lubrication

3.2 Routine check-ups

3.3 Cleaning

3.4 Storage

4. Technical specification

4.1 Dimensions

4.2 Labels

5. EU-Declaration of Conformity

6. Contact details

FOREWORD

Welcome as the owner of a product manufactured for Tarsus by Sjöbloms Sjukvårdsutrustning AB. We hope that you will be pleased and satisfied with our quality product.

If a problem should arise with this or any of our products, you are always welcome to contact us, and our service department will assist you with the expertise needed to solve the problem.

The user manual contains all the information needed to use, maintain and service the product. The document is continuously updated as the products are modified. Therefore we reserve against possible changes in the user manual.

WARRANTY

Sjöbloms provides a two-year warranty for material, design and manufacturing defects, provided the products are not subjected to tampering or excessive wear. For the guarantee to apply, the product must have been used in accordance with the user manual, and the maintenance performed as specified.

IMPORTANT INFORMATION

This product is intended for professional use only. The user manual should be read carefully before using the product, to avoid unnecessary risks for the patient and operator.

The user manual contains warnings indicating risk of personal injury, which should be noted. Improper handling increases the risk.

A warning triangle is shown where particularly important information can be found.



Warning triangle

1. DESCRIPTION OF FUNCTIONS

The Airport emergency transport trolley is intended for professional use only, for transportation of a patient between ambulance and airplane.

It is designed to be a functional and user-friendly tool for the operator and patient.

Note that the product must only be handled by staff who have read the user-manual, so that it is handled correctly.

General description

Airport emergency transport trolley, hydraulic height adjustment, coated all-welded steel chassis.

200 mm central locking castors, two of which has a directional locking function available.

2. OPERATION

2.1 Height adjustment

Height is raised by pressing the pedal (A). It is lowered by turning the knob (B) counter-clockwise.



Note!

Risk of pinching. Be cautious when adjusting the height.

2.2 The castors

The trolley has 200 mm central locking castors.

The castors are locked when the red pedal is pressed down.

The directional-locking function is available as option.



2.3 The Slide


The red stops prevents the slideboard to move more than intended.



Note!

Risk of pinching.

To free the slide, pull out the lock (marked) and pull out the slide in the arrow's direction.



Emergency release of the stretcher

If the normal release of the stretcher fails for some reason, the stretcher may be released manually. Press both the buttons inwards towards the centre of the trolley.



2.4 Against the airplane

Make sure that the top of the trolley is in its start position, marked as shown arrow against arrow. If it's not, release the top by turning the handle and adjust into position.

Push the trolley towards the airplane gently. Remember that the trolley is heavy. Use low speed and stop some inches before the airplane to avoid damage. Lock the castors by pressing the red part of the pedal (A).




NOTE! The castors must be locked before transfer of patient.



Adjust the height by pressing B until it fits the airplane height. Release the slide and pull out. Make sure the slide height is slightly higher than the airplane attachment spot.

Match the slot (C) of the moving attachment plate (D) to the matching airplane ridge. Lower the trolley gently so that the slot is locked in the ridge.



NOTE! Don't lower the trolley more than necessary to lock the slot/ridge, otherwise the slide may brake.

When the attachment plate of the slide is in place, the top must be locked by turning the handle downwards.




NOTE!

The patient must not be transferred to the trolley unless the slide is in right position and the top is locked.



When the patient is moved from the airplane to the trolley, the slide should be pushed-in until locked by the fast-lock.



NOTE!

Lower the trolley to a low position (minimum below the handles) before the trolley is moved with a patient.

If lowered to its absolute lowest position, a slight tremble may occur that could be troubling for the patient.



3. MAINTENANCE AND CLEANING

3.1 Lubrication

During maintenance and repair of hydraulic, the red safety locks must be activated (2 pcs).


The trolley must be raised to mid-height for the locks to be activated.



Lubricate all joints. Use a thin lubrication oil, 1-3 drops, 1-2 times per month.




Grease the excenter contact space once a month.



Check oil level.
Volume is 1.1 litre.



Inspect tubes and pump for
leakage or cracks.



Roller surface. Lubricate all rolls with a
thin oil.



3.2 Routine controls

Routine controls and lubrication is recommended for every 6 months or when needed to maintain a fully functioning trolley.



IMPORTANT!

Failure to perform the following routine can cause the patient and/or operator to be exposed to unnecessary risks and the trolley to have a reduced life expectancy. Please check all functions and parts on a regular basis.

Check the following:

- Any unusual sounds.
- Hydraulic pump function
- Castors and castor locking mechanism.
- That joints and fittings are lubricated according to the lubrication chart.

3.3 Cleaning

Can be cleaned and disinfected.

3.4 Storage

The trolley is manufactured and designed to withstand moist and snow. To maintain functions of castors, locking mechanisms and the slide however, we recommend that the trolley is stored protected, however not necessarily in-doors.

4. TECHNICAL SPECIFICATIONS

4.1 Dimensions

Length

2400 mm

With slide pulled-out: 3800 mm

Width

900 mm.

Height

Max 1700 mm

Min 550 mm


Weight

175 kg

Max load

300 kg

4.5 Labels



The CE label contains the date, order number, serial number and our contact information.



Max load 300 kg.



Risk of pinching



5. DECLARATION OF CONFORMITY

Type of equipment: Airport emergency transport trolley

Art no: B6-030-4000

Product reg no: 60211

Manufacturer

Sjöbloms Sjukvårdsutrustning AB

Myrängsvägen 8

S-891 50 ÖRNSKÖLDSVIK

SWEDEN


Sjöbloms Sjukvårdsutrustning AB ensures that the above products are in accordance with the EU Council MDD 93/42/EEC, as well as the Swedish legislation in force (1993:584) on medical devices and the regulations issued by Medical Products Agency (LVFS 2003:11).

Örnsköldsvik 2015-08-28

David Sjöblom

CEO

6. CONTACT INFO



Sjöbloms Sjukvårdsutrustning AB
Myrängsvägen 8
891 50 Örnsköldsvik
Sweden

Tel: +46 660 105 80

Fax: + 46 660 180 00

E-mail: info@sjobloms.com

Hemsida: www.sjobloms.com