

## Undersökningslampor

### Halux® 50 P SX

090-5000 Halux 50 P SX, 2-steps strömbrytare.

#### Tillbehör

##### Art.nr.

090-5001 Väggfäste för Halux 50 P SX  
090-5002 Halogenlampa Halux  
090-5003 Golvstativ för Halux 50 P SX  
090-5004 Fäste för Halux 50 P SX, för plac. på gynstol  
090-5005 IVA-fäste för Halux 50 P SX

#### Ljusstyrka:

Ljusstyrka avstånd 0,5 m, ca 48 000 lux

#### Vikt:

Ca 1,4 kg

#### Glödlampa

Effekt lampa 50 W



090-5000 och 090-5004  
på vår gynstol 015-2031



090-5003 Golvstativ



090-5001 Väggfäste



090-5005 IVA-fäste

# Dmed® halux 50 P SX

## TECHNICAL DATA

## EXAMINATION LIGHTS

FITTED WITH	1 x halogen lamp QR-CBC51 50 W - D68727000
LUMINOUS POWER	48 000 Lux / 0,5 m
WORK EQUIPMENT	electronic transformer
CONNECTED LOAD	230 V; 50 Hz
MAINS LEAD	approx. 3 m Euro
LUMINAIRE BODY	polycarbonate (PC)
WEIGHT (NET)	approx. 1,4 kg
POWER CONSUMPTION	approx. 50 W
USAGE	rocker switch (I/O), 2-stages rocker switch (60/100 %)
TECHNOLOGY	2-pole primary
BALANCE OF ARTICULATION JOINTS	friction
CLASS OF PROTECTION	II
ACCESSORY	colour correction filter (ca. 4100 K)
NORMS	EN 60601-1, EN 60598-2-25
FASTENING	pin (16 mm W), rail, roller stand, table, wall
ORDER NUMBER	D13371000 - pure white



## LIGHT DISTRIBUTION

## TECHNICAL DRAWING

