

PERI



## Tested cleaning products

- Always use the disinfectant according to the instructions of the manufacture.
- The disinfectants should not be changed.
- Only use one disinfectant and rinse with clear water; never use a mixture of various disinfectants.
- As we have no influence on the production and raw materials of the disinfectants, we can give no general guarantee of resistance.
- We recommend testing a disinfectant on an invisible point of the covering before any regular application of these disinfectants.
- Products marked with a \* are quick-acting disinfectants.

MANUFACTURER	DISINFECTANT
Antiseptica Dr. Hans-Joachim Molitor	Antiseptica Kombi Flächen Des.
	Biguacid Plus
	Biguacid-S
	Descogen® Liquid
	Descocid®-N
B. Braun Melsungen AG	Meliseptol® rapid *
	Melsept® SF
	Hexaquart® S
	Meliseptol® foam pure + fresh
Bode Chemie GmbH	Bacillo® AF *
	Bacillo® AF-Tissues *
	Bacillo® 30 Foam
	Bacillo® 30 Tissues *
	Kohrsolin® extra Tissues
	Mikrobac® forte
	Mikrobac® Tissues *
CosiMed GmbH	CosiMed Schnelldesinfektion *
	CosiMed Desinfektionstücher
DESOMED	DESOMED® RAPID UNIVERSAL *
DR. SCHUMACHER GMBH	BIGUANID FLÄCHE N
	CLEANISEPT
	DESCOSEPT SENSITIVE WIPES *
	DESCOSEPT AF *
	OPTISAL® N
	OPTISEPT®
DR. Weigert	Neoform MED rapid *

MANUFACTURER	DISINFECTANT
Ecolab Deutschland GmbH	Incidin® Active
	Incidin® Foam *
	Incidin™ OxyWipe
	Incidin® PLUS
	Incidin® Liquid Spray *
	Sani-Cloth® Active – Tücher *
GAMA Healthcare Ltd.	Clinell Alcohol Wipes
	Clinell Universal Wipes
JOYA International	JOYA®
Lysoform	Aldasan® 2000
	Clorina®
	Fugaten® *
	Hospisept® - Tuch *
	Lysoformin®
	Lysoformin® spezial
	Trichlorol®
OPTIMUM11 GMBH	BESTWIPES *
Peppler GmbH	Micro Clean *
	Micro Clean Foam *
	Immix Clean *
PINO GmbH	Septapin® Wipes *
Reckitt Benckiser	SAGROTAN® Hygiene-Spray
SARAYA EUROPE	Alpet D2 *
Schülke & Mayr GmbH	Mikrozyd® AF liquid *
	Mikrozyd® AF wipes *
	Mikrozyd® sensitive liquid *
	Mikrozyd® sensitive wipes *
	Perform®
	Terralin® liquid *
	Terralin® protect
	TPH protect
Tana Chemie GmbH	APESIN rapid